

[Skin camera]

# Charm Lite Fixed Lens

50×PL

To examine skin texture or pigmentation



200×

To examine the condition of the hair or scalp



*An economical model that is easy to operate*  
*A hand-held model that achieves maximum performance with simple operations.*



You do not need to focus or adjust magnification.  
 You can operate this scope with one simple switch.  
 Connect the USB cap, and transfer skin images acquired by the Charm Lite into the personal computer.

PC

USB-CAP




Charm Lite

■ Specifications

Camera (built-in)	Image pickup device:	1/3-inch CCD color image sensor
	Effective pixels:	510(H)×492(V) equivalent to 250,000-pixel
	Scanning method:	2:1 interlace
	White balance:	Fixed
Light source	AGC:	ON
	LED:	White LED
	Light control:	Fixed
	Lens:	Fixed, Charm Lite III 50PL: 50× Charm Lite III 200PL: 200×
	Cooling:	Air-cooling method
Dimensions:	50(W)×L(D)×165(H) (50PL : L=70/200 : L=80)	
Weight:	About 140 g	

Stand	Image output:	VBS 1Vp-p, conform to NTSC: 1 system
	Image input:	Charm Lite III: 1 system VBS 1Vp-p, conform to NTSC: 1 system
	Dimensions:	70(W)×150(D)×90(H)
	Weight:	About 200g
Accessories	Input power supply:	DC 12 V 800 mA
	Working temperature:	5°C - 35°C
	Working humidity:	35% - 85% (no condensation)
	AC adaptor:	AC adaptor: 1 unit
	Curl cord:	Curl cord: 1
	Video cord:	Video cord: 1
	Operation manual:	Operation manual: 1

■ Dimensions (mm)



\*The device does not have static image or frame memory functions.